

**WEEKLY**

CPAP SOAP PACKETS

Mix 5 ml soap packet with 1 gallon of water to clean CPAP accessories. Rinse products thoroughly. Allow the product to air dry. Follow CPAP manufacturers guidelines for cleaning.



CPAP SOAP BUBBLE PADS

Use pads on CPAP accessories and discard the pad after cleaning. Allow CPAP accessories to air dry. Use a clean cloth to remove any excess bubbles. If irritation occurs, stop use and seek medical advice. Follow CPAP manufacturer's guidelines for cleaning. Store pads in a cool dry location.

**DAILY**